

Filter Wheel Testing Report

Project No.:	17125
End User:	University of California, Santa Cruz
End User Country:	United States

Manufacturer	Sciencetech Inc.
Filter Wheel Model:	SFW-6-1
System Voltage Specification:	24VDC

Performed By:	JF
Date [YYYY-MM-DD]:	2024-09-11
Reviewed By:	RP
Date of Issue [YYYY-MM-DD]:	2024-09-12

Order Sort Filters and Their Location on the Wheel with Respect to Home Position	
Filter at Position 1:	No Filter Located at this position
Filter at Position 2:	Long Pass Filter – Cut-On Wavelength: 295nm
Filter at Position 3:	Long Pass Filter – Cut-On Wavelength: 500nm
Filter at Position 4:	Long Pass Filter – Cut-On Wavelength: 900nm
Filter at Position 5:	Long Pass Filter – Cut-On Wavelength: 1500nm
Filter at Position 6:	No Filter Located at this position

Notes:

1. It is important to mount the filter wheel in the beam path with the correct orientation, as indicated by the transmission arrow placed on its housing. Otherwise, the filters may be permanently damaged by light passing through in the wrong direction.
2. A given filter position number indicates which filter will be located in the beam path when that position is selected.